



February 2020



Helpful Tips for Using Your CapTel Phone

Tip: Take Your Time Dialing

With CapTel, you can dial the phone number of the person you are calling BEFORE you pick up the handset. This gives you extra time to make sure the phone number is right before the call goes through.

Simply leave the CapTel handset down (hung up) and press the phone number buttons on the dial keypad. Don't forget to press a "1" or the area code if needed. As you dial, you can see the phone number on the CapTel display screen. If you make a mistake, erase the number by pressing the **DOWN ARROW** button on the CapTel 840 / 840i / 880i or touching Backspace (**X**) on the CapTel 2400i.



When you are ready to place your call, pick up the CapTel handset. Your phone will dial the number you entered and the call will go through.



No More Repeats

As Eleanor's hearing loss became more noticeable over time, she grew tired of asking people to repeat themselves over the phone. "I was forever saying, I'm sorry what did you say? Please repeat that? That kind of thing," she explains. A friend used CapTel and recommended that Eleanor try it herself. "It's been the answer," she exclaims. Eleanor especially relies on CapTel to keep up with family activities. "I wouldn't be without it," she says.

See Eleanor's story »



I love my CapTel phone and I am thankful for it every day. I would not be able to use the telephone at all without it! - CapTel user in Connecticut I love this phone, I really do! - CapTel user in Texas I love my CapTel phone. It makes it easy to use the phone.

- CapTel user in Florida

'Round the Clock Help – Just one button away! Help is available 24 hours a day, 7 days a week! (*Excluding major holidays.*)

To contact:

6

0

- Just press the blue button on your CapTel phone!
- 1-888-269-7477 1-866-670-9134 (Spanish)
- Online Chat »

CapTel@CapTel.com »

Connect with us online!



Thank you to everyone who follows CapTel on Facebook. With more than 38,000 followers, the online community of CapTel users is another great resource.

Proud Sponsor of



Have a CapTel story to share? Send it to Newsletter@CapTel.com. We'd love to hear from you!

CapTel.com » | CapTel Support » | CapTel News & Blog »

CapTel Inc. 450 Science Drive, Madison, WI 53711 © 2020 Ultratec, Inc.

THE CAPTIONS TURNED ON. IP Captioned Telephone Service may use a live operator. The operator generates captions of what the other party to the call says. These captions are then sent to your phone. There is a cost for each minute of captions generated, paid from a federally administered fund. No cost is passed on to the CapTel user for using the service.

CapTel captioning service is intended exclusively for individuals with hearing loss. CapTel® is a registered trademark of Ultratec, Inc. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Ultratec, Inc. is under license. Other trademarks and trade names are those of their respective owners. (v10.1 10-19)