




Updating Your *CapTel* Phone

Occasionally, changes are made to the *CapTel* 800 software. Your phone can check whether new software is available and then perform an update as needed over the telephone line. This process ensures you will always have the latest technology for as long as you own your *CapTel* 800.

To update your *CapTel* 800 software:

-  1. With the handset hung up, press the **YES** button to see Options.
- 
 2. Press the **DOWN** arrow button repeatedly until **Settings** is highlighted. Press the **YES** button to select.
- 
 3. Press the **DOWN** arrow button repeatedly until **Phone Settings** is highlighted. Press the **YES** button to select.
- 
 4. Press the **DOWN** arrow button repeatedly until **Update Phone** is highlighted. Press the **YES** button to select.





5. Press the **YES** button to begin updating your phone, or press the **NO** button to exit.



NOTE: *You can set the update to begin automatically at a more convenient time. Press the UP arrow button to set up a timed delay. Select when you want the update to begin (in 1 hour, in 2 hours, etc.) using the UP and DOWN arrow buttons. Once a delay has been set, you can continue using your CapTel phone to make/answer calls. This will not disrupt the planned update process.*

6. The *CapTel* phone checks to make sure that your software is up to date. If a new software version is available, your *CapTel* phone automatically updates to the most current updates to the most current version, then resets itself so you are ready to make calls.

